



M2 Automation
Eresburgstr. 22 ♦ 12103 Berlin

Pump Unit 002

PUU02

A universal liquid pump unit with adjustable pump speed.



Indicator

Pump Unit PUU02 A/B/C/D

- Max flow rate
 - A 0.05 l/min
 - B 0.1 l/min
 - C 0.3 l/min
 - D 0.6 l/min
- Chemically inert (see notes)
- Adjustable pump speed (option)
- Pulsation Damper (option)
- 24V DC motor (12V Version on request)
- Dimensions 165x110x55mm (LxWxH)
- 24V DC power adapter



Power in

Luer-Lock port in

Luer-Lock port out

Notes:

Liquid will have contact with the listed materials below:

EPDM = Ethylene Propylene Diene Monomer; PP = Polypropylene; PTFE = Poly(tetrafluoroethylene)

M2 Automation
Martin Mueller
Eresburgstr. 22
12103 Berlin
Deutschland

Mobil
Telefon
Fax
E-Mail
Internet

+49 (0) 173.62 05 717
+49 (0) 30.53 67 11- 44
+49 (0) 30.53 67 11- 45
info@m2-automation.de
www.m2-automation.de